



TREVALLI

TREVALLI SURGICAL

Quality is our Priority

ORTHOPEDIC IMPLANTS

INSTRUMENTS

2025



3 Oldham Place, Pietermaritzburg



083 555 8643



trevalli@mweb.co.za

Table Of Content

LOCKING PLATES

Page No.1-11

LOCKING SCREW

Page No.12-13

BONE PLATE

Page No.14-22

BONE SCREW

Page No.23-27

INTERLOCKING NAIL SYS

Page No.28-31

INTERLOCKING BOLT SYS

Page No.32-33

PROSTHESIS SYSTEM

Page No.34

SPINE SYSTEM

Page No.35-39




TREVALLi

LOCKING PLATES



Small Compression Locking Plate 3.5mm

Screw	3.5mm
Holes	3 to 16
Material	Stainless Steel 316L,316LVM Titanium
Shaft Length	10mm to 40mm (diff 2.0mm)




LCP Small T Plate (Oblique) 3.5mm - Left & Right

Screw	3.5mm
Holes	2 to 10
Material	Stainless Steel 316L,316LVM Titanium
Shaft Length	10mm to 40mm (diff 2.0mm)



LCP Proximal Humerus Plate 3.5mm (Type 1)

Screw	3.5mm & 4.0mm
Holes	4 to 12
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (diff 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (diff 2.0mm)



LCP Small T Plate (Oblique) 3.5mm - Right Angle

Screw	3.5mm
Holes	2 to 10
Material	Stainless Steel 316L,316LVM Titanium
Shaft Length	10mm to 40mm (diff 2.0mm)



TREVALLi

LOCKING PLATES



LCP Proximal Humerus Plate 3.5mm (Type 2)

Screw	3.5mm & 4.0mm
Holes	3 to 10
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (diff 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (diff 2.0mm)

LCP Proximal Humerus Plate 3.5mm (Type 3)

Screw	3.5mm & 4.0mm
Holes	3 to 8
Material	Stainless Steel 316L, Titanium
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (diff 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (diff 2.0mm)



LCP Cloverleaf Plate 3.5mm

Screw	3.5mm
Holes	3 to 13
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (diff 2.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (diff 2.0mm)

LCP Cloverleaf Plate 3.5mm

Screw	3.5mm
Holes	3 to 18
Material	Stainless Steel 316L,316LVM Titanium
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (diff 2.0mm)



TREVALLi

LOCKING PLATES



LCP Reconstruction Plate 5.0mm (Round hole)

Screw	5.0mm
Holes	3 to 18
Material	Stainless Steel 316L, 316LVM Titanium
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)

LCP Reconstruction Plate 3.5mm (Dynamic Hole)

Screw	3.5mm
Holes	3 to 16
Material	Stainless Steel 316L, 316LVM Titanium
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (diff 2.0mm)



Distal Tibia Locking Plate Left & Right

Screw	4.0mm & 5.0mm
Holes	4 to 12
Material	Stainless Steel 316L, 316LVM Titanium
Proximal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)



LCP Reconstruction Plate 5.0mm (Dynamic Hole)

Screw	5.0mm
Holes	3 to 16
Material	Stainless Steel 316L, 316LVM Titanium
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)



TREVALLi

LOCKING PLATES



Medial Distal Tibia Locking Plate 3.5mm Left & Right

Screw	3.5mm & 4.0mm
Holes	5 to 13
Material	Stainless Steel 316L, 316LVM Titanium
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (diff 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (diff 2.0mm)



Distal Tibia Anatomical Locking Plate 3.5mm Left & Right

Screw	3.5mm & 4.0mm
Holes	4 to 14
Material	Stainless Steel 316L, 316LVM Titanium
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (diff 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (diff 2.0mm)



Narrow Locking Compression Plate 5.0mm

Screw	5.0mm
Holes	4 to 18
Material	Stainless Steel 316L, 316LVM Titanium
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)



Broad Locking Compression Plate 4.8mm

Screw	5.0mm
Holes	4 to 20
Material	Stainless Steel 316L, 316LVM Titanium
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)



TREVALLI

LOCKING PLATES




T Plate Locking 5.0mm

Screw	5.0mm
Holes	4 to 12
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)



L Butterss Locking Plate Left & Right

Screw	5.0mm
Holes	4 to 12
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)



T Butterss Locking Plate 5.0mm

Screw	5.0mm
Holes	4 to 12
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)



Distal Humerus Locking Plate Left & Right

Screw	3.5mm
Holes	3 to 15
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)



TREVALLi

LOCKING PLATES



Distal Humerus with Support Locking Plate Left & Right

Screw	3.5mm
Holes	3 to 15
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)



Extra-Articular Distal Humerus Locking Plate Left & Right

Screw	3.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)



Lateral Tibial Locking Plate 5.0mm/6.5mm - Left & Right

Screw	5.0mm
Holes	3 to 13
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Cancellous Locking Screw 5.0mm
Screw Length	30mm to 80mm (diff 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)



Lateral Tibial Locking Plate 5.0mm/6.5mm - Left & Right

Screw	5.0mm & 6.5mm
Holes	3 to 12
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Locking Screw 6.5mm
Screw Length	30mm to 95mm (diff 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)



TREVALLi

LOCKING PLATES



Distal Femur Locking Plate 5.0m/6.5mm (Left & Right)

Screw	5.0mm & 6.5mm
Holes	3 to 13
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Cancellous Locking Screw 6.5mm
Screw Length	30mm to 95mm (diff 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)



Distal Femur Locking Plate Left & Right

Screw	5.0mm
Holes	3 to 13
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Cancellous Locking Screw 5.0mm
Screw Length	30mm to 80mm (diff 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)



Condylar Butters Locking Plate Left & Right

Screw	5.0mm & 6.5mm
Holes	3 to 13
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Cancellous Locking Screw 6.5mm
Screw Length	30mm to 95mm (diff 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)



DHS Plate locking long Barral^o130 /135

Screw	5.0mm & 6.5mm
Holes	2 to 12
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)



TREVALLi

LOCKING PLATES



D.C.S Plate Locking 90

Screw	5.0mm & 6.5mm
Holes	3 to 16
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)
Screw 6.5mm	30mm to 95mm (diff 5.0mm)

DHS Plate locking Short Barra 130 /135

Screw	5.0mm
Holes	2 to 12
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Locking Screw 5.0mm
Screw Length	10mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)



Proximal Femur Locking Plate

Screw	5.0mm & 6.5mm
Holes	5 to 13
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 5.0mm
Screw Length	10mm to 0mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)
Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 100mm (diff 5.0mm)

Distal Femur Locking Plate 5.0mm/6.5mm(Left & Right)

Screw	3.5mm
Holes	4 to 8
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 3.5mm
Screw Length	12mm to 40mm (diff 2.0mm)



TREVALLI

LOCKING PLATES



Proximal Tibia Locking Plate Left & Right

Screw	5.0mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 3.5mm
Screw Length	12mm to 40mm (diff 2.0mm)



LCP Tomfix Mid High T Plate

Screw	5.0mm
Holes	03, 04
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 5.0mm
Screw Length	18mm to 50mm (diff 2.0mm)
Screw Length	55mm to 90mm (diff 5.0mm)



Distal Tibia with Tab Locking Plate Left & Right

Screw	3.5mm & 4.0mm
Holes	5 to 13
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (diff 2.0mm)
Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (diff 5.0mm)




Calcaneus Locking Plate Left & Right

Screw	3.5mm
Holes	70mm, 75mm
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 3.5mm
Screw Length	12mm to 40mm (diff 2.0mm)



TREVALLi

LOCKING PLATES



2.4mm Distal Radius Locking Plate Straight

Screw	2.4mm
Holes	5 to 13
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 2.4mm
Screw Length	06mm to 30mm (diff 2.0mm)



2.4mm Distal Radius Locking Plate (Head 2 Hole) Left & Right

Screw	2.4mm
Holes	2+3, 2+4
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 2.4mm
Screw Length	06mm to 30mm (diff 2.0mm)



2.4mm Distal Radius Locking Plate (Head 3 Hole) Left & Right

Screw	2.4mm
Holes	3+3, 3+4
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 2.4mm
Screw Length	06mm to 30mm (diff 2.0mm)



2.4mm L Distal Radius Oblique Locking Plate Left & Right

Screw	2.4mm
Holes	3+3, 3+4
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 2.4mm
Screw Length	06mm to 30mm (diff 2.0mm)



TREVALLi

LOCKING PLATES



2.4mm Distal Radius Locking Plate (Head 5 Hole) Left & Right

Screw	2.4mm
Holes	5+3, 5+5
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 2.4mm
Screw Length	06mm to 30mm (diff 2.0mm)



2.4mm Distal Radius Locking Plate (Head 3 Hole)

Screw	2.4mm
Holes	3+3, 3+4
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 2.4mm
Screw Length	06mm to 30mm (diff 2.0mm)



2.4mm Buttress Locking Plate (Head 5 Hole) Left & Right

Screw	2.4mm
Holes	5+3
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 2.4mm
Screw Length	06mm to 30mm (diff 2.0mm)



Two-Column Distal Radius Plate Left & Right

Screw	2.4mm
Holes	2 to 4
Material	Stainless Steel 316L,316LVM Titanium
Screw	Locking Screw 2.4mm
Screw Length	06mm to 30mm (diff 2.0mm)



TREVALLi

LOCKING SCREW



LCP Screw DIA 2.4mm

Thread Diameter	2.4mm
Use Drill Diameter	2.0mm
Material	Stainless Steel 316L,316LVM Titanium
Drive	Hex
Screw Length	06mm to 30mm (diff 2.0mm)
Size	2.4mm




LCP Screw DIA 2.7mm

Thread Diameter	2.7mm
Use Drill Diameter	2.0mm
Material	Stainless Steel 316L,316LVM Titanium
Drive	Hex
Screw Length	06mm to 30mm (diff 2.0mm)
Size	2.7mm



LCP Screw DIA 3.5mm

Thread Diameter	4.0mm
Use Drill Diameter	3.2mm
Material	Stainless Steel 316L,316LVM Titanium
Drive	Hex
Screw Length	25mm to 70mm (diff 2.0mm)
Size	4.0mm




LCP Screw DIA 3.5mm

Thread Diameter	3.5mm
Use Drill Diameter	2.7mm
Material	Stainless Steel 316L,316LVM Titanium
Drive	Hex
Screw Length	10mm to 50mm (diff 2.0mm)
Size	3.5mm



LCP Screw DIA 5.0mm

Thread Diameter	5.0mm
Use Drill Diameter	4.0mm
Material	Stainless Steel 316L,316LVM Titanium
Drive	Hex
Screw Length	16mm to 50mm (diff 2.0mm)
Screw Length	55mm to 95mm (diff 5.0mm)
Size	5.0mm



LCP Screw Cancellous DIA 3.5mm

Thread Diameter	5.0mm
Use Drill Diameter	4.0mm
Material	Stainless Steel 316L,316LVM Titanium
Drive	Hex
Screw Length	30mm to 50mm (diff 2.0mm)
Size	5.0mm



TREVALLI

LOCKING SCREW



LCP Screw DIA 6.5mm

Thread Diameter	6.5mm
Use Drill Diameter	4.5mm
Material	Stainless Steel 316L,316LVM Titanium
Drive	Hex
Screw Length	35mm to 95mm (diff 2.0mm)
Size	6.5mm

TREVALLI



TREVALLI

BONE PLATE



Dynamic Compression Plate for 2.7mm Screw

Screw	2.7mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 2.7mm
Screw Length	06mm to 30mm (diff 2.0mm)



Reconstruction Plate for 2.7mm Screw Round Hole

Screw	2.7mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 2.7mm
Screw Length	06mm to 30mm (diff 2.0mm)



Reconstruction Plate for 3.5mm Screw Round Hole

Screw	3.5mm
Holes	4 to 18
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)



One Third Tubular Plate For 3.5mm Screw

Screw	3.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)



Reconstruction Curved Plate for 3.5mm Screw

Screw	3.5mm
Holes	4 to 16
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)



Reconstruction Plate for 3.5mm Screw DC Hole

Screw	3.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)



TREVALLi

BONE PLATE



LC DCP Plate for 3.5mm Screw

Screw	3.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)



Small DCP Plate for 3.5mm Screw

Screw	3.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)



Small T Oblique Plate for 3.5mm Screw

Screw	3.5mm
Holes	3 to 10
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)



Small T Plate for 3.5mm Screw

Screw	3.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)



TREVALLi

BONE PLATE



Acro Plate for 3.5mm Screw Left & Right

Screw	3.5mm
Holes	04,05,06
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)



Semi Tubular Plate for 4.5mm

Screw	4.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)



Reconstruction Plate for 4.5mm Screw Round Hole

Screw	4.5mm
Holes	4 to 18
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)



Reconstruction Plate for 4.5mm Screw DC Hole

Screw	4.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)



TREVALLi

BONE PLATE



LC DCP Plate Narrow for 4.5mm Screw

Screw	4.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)



Narrow DCP Plate for 4.5mm Screw

Screw	4.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)



Broad DCP Plate for 4.5mm

Screw	4.5mm
Holes	4 to 18
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)



LC DCP Plate Broad for 4.5mm Screw

Screw	4.5mm
Holes	6 to 14
Material	Stainless Steel 316L,316LVM Titanium
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)

BONE PLATE



T Plate for 4.5mm Screw

Screw	6.5mm & 4.5mm
Holes	3 to 12
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (diff 2.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)



Cloverleaf Plate

Screw	
Holes	3 to 12
Material	Stainless Steel 316L,316LVM Titanium
Screw	
Screw Length	



L Buttress Plate for 4.5mm Screw Left & Right

Screw	6.5mm & 4.5mm
Holes	3 to 12
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (diff 2.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)



T Buttress Plate for 4.5mm Screw

Screw	6.5mm & 4.5mm
Holes	3 to 12
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (diff 2.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)



TREVALLI

BONE PLATE



L Buttress Plate for 4.5mm Screw Left & Right

Screw	6.5mm & 4.5mm
Holes	3 to 12
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (diff 2.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)

Condylar Buttress Plate for 4.5mm Screw Left & Right

Screw	6.5mm & 4.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (diff 2.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)



DHS Plate (Round Holes) Short Barrel

Screw	12.5mm & 4.5mm
Holes	4 to 18
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Leg Screw 12.5mm
Screw Length	50mm to 115mm (diff 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)
Degree	120°, 125°, 130°, 135°, 140°, 145°

DHS Plate (DCP Holes) Short Barrel

Screw	12.5mm & 4.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Leg Screw 12.5mm
Screw Length	50mm to 115mm (diff 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)
Degree	120°, 125°, 130°, 135°, 140°, 145°



TREVALLI

BONE PLATE



DHS Plate (DCP Holes) Long Barrel

Screw	12.5mm & 4.5mm
Holes	2 to 18
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Leg Screw 12.5mm
Screw Length	50mm to 115mm (diff 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)
Degree	120°, 125°, 130°, 135°, 140°

DHS Plate (DCP Holes) Long Barrel

Screw	12.5mm & 4.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Leg Screw 12.5mm
Screw Length	50mm to 115mm (diff 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)
Degree	120°, 125°, 130°, 135°, 140°, 145°



DHS Plate (Round Holes) Pediatrics

Screw	12.5mm & 4.5mm
Holes	4 to 10
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (diff 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)
Degree	125°, 130°, 135°

950DCS Plate (DCP Holes)

Screw	12.5mm & 4.5mm
Holes	6 to 20
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (diff 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)



TREVALLi

BONE PLATE



Y Reconstruction Plate for 3.5mm Screw

Screw	3.5mm
Holes	1 to 10
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (diff 2.0mm)
Distal Screw	Cancellous Screw 4.0mm
Screw Length	10mm to 30mm (diff 2.0mm)
Distal Screw	Cancellous Screw 4.0mm
Screw Length	35mm to 70mm (diff 5.0mm)



DCS Plate (Round hole) 95

Screw	12.5mm & 4.5mm
Holes	6 to 20
Material	Stainless Steel 316L,316LVM Titanium
Proximal Screw	Leg Screw 12.5mm
Screw Length	50mm to 115mm (diff 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)
Distal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (diff 5.0mm)



Cobra Head Plate for 4.5mm Screw Left & Right

Screw	4.5mm & 6.5mm
Holes	4 to 14
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Leg Screw 4.5mm
Screw Length	16mm to 70mm (diff 2.0mm)
Distal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (diff 5.0mm)



Calcaneal Plate 3.5mm

Screw	4.0mm & 3.5mm
Holes	Small Medium Large
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Cortical Screw 3..5mm
Screw Length	10mm to 60mm (diff 2.0mm)
Distal Screw	Cancellous Screw 4.0mm
Screw Length	10mm to 30mm (diff 2.0mm)



TREVALLI

BONE PLATE



Trocentric Plate For DHS

Screw	4.5mm & 6.5mm
Holes	4 Holes
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Cortical Screw 4..5mm
Screw Length	16mm to 70mm (diff 2.0mm)
Distal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (diff 5.0mm)



Eight Plate

Screw	4.5mm
Holes	4 Holes
Material	Stainless Steel 316L,316LVM Titanium
Distal Screw	Cortical Screw 4..5mm



TREVALLi

BONE SCREW



Cortical Screw 2.7mm

Material	Stainless Steel 316L,316LVM Titanium
Size	2.7mm
Length	06mm to 30mm (diff 2.0mm)



Cortical Screw (Self Tapping) 1.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	1.5mm
Length	06mm to 30mm (diff 2.0mm)



Cortical Screw (14 TPI) 3.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	3.5mm
Length	10mm to 60mm (diff 2.0mm)



Cortical Screw 2.0mm

Material	Stainless Steel 316L,316LVM Titanium
Size	2.0mm
Length	06mm to 30mm (diff 2.0mm)



Cortical Screw (14 TPI) Self Taping 3.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	3.5mm
Length	10mm to 60mm (diff 2.0mm)



Cortical Screw 2.4mm

Material	Stainless Steel 316L,316LVM Titanium
Size	2.4mm
Length	06mm to 30mm (diff 2.0mm)



TREVALLi

BONE SCREW



Cortical Screw (20 TPI) 3.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	3.5mm
Length	10mm to 60mm (diff 2.0mm)



Cortical Screw (20 TPI) Self Tapping 3.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	3.5mm
Length	10mm to 60mm (diff 2.0mm)



Cancellous Screw (Full Thread) 4.0mm

Material	Stainless Steel 316L,316LVM Titanium
Size	4.0mm
Length	10mm to 30mm (diff 2.0mm)
Length	35mm to 70mm (diff 5.0mm)



Cancellous Screw (Short Thread) 4.0mm

Material	Stainless Steel 316L,316LVM Titanium
Size	4.0mm
Length	10mm to 30mm (diff 2.0mm)
Length	35mm to 70mm (diff 5.0mm)



Cannulated Screw (Short Thread) 4.0mm

Material	Stainless Steel 316L,316LVM Titanium
Size	4.0mm
Length	10mm to 30mm (diff 2.0mm)
Length	35mm to 70mm (diff 5.0mm)



Cannulated Screw (Full Thread) 4.0mm

Material	Stainless Steel 316L,316LVM Titanium
Size	4.0mm
Length	10mm to 30mm (diff 2.0mm)
Length	35mm to 70mm (diff 5.0mm)



TREVALLi

BONE SCREW



Cortical Screw 4.5mm DIA

Material	Stainless Steel 316L,316LVM Titanium
Size	4.5mm
Length	16mm to 70mm (diff 2.0mm)



Cortical Screw Self Tapping 4.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	4.5mm
Length	16mm to 70mm (diff 2.0mm)



Cancellous Screw (16th) 6.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	6.5mm
Length	20mm to 110mm (diff 5.0mm)



Malleolar Screw 4.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	4.5mm
Length	20mm to 100mm (diff 5.0mm)



Cancellous Screw (Full TH) 6.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	6.5mm
Length	25mm to 110mm (diff 5.0mm)



Cancellous Screw (32th) 6.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	6.5mm
Length	20mm to 110mm (diff 5.0mm)



TREVALLi

BONE SCREW



Cannulated Screw (32th) 6.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	6.5mm
Length	25mm to 110mm (diff 5.0mm)




Cannulated Screw (16th) 6.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	6.5mm
Length	20mm to 110mm (diff 5.0mm)



DHS Lag Screw

Material	Stainless Steel 316L,316LVM Titanium
Size	12.5mm
Length	50mm to 110mm (diff 5.0mm)



Cannulated Screw (Full Thread) 6.5mm

Material	Stainless Steel 316L,316LVM Titanium
Size	6.5mm
Length	20mm to 100mm (diff 5.0mm)



DHS Lag Screw (Hex)

Material	Stainless Steel 316L,316LVM Titanium
Size	12.5mm
Length	50mm to 110mm (diff 5.0mm)



Harbat Screw

Material	Stainless Steel 316L,316LVM Titanium
Size	4.0mm
Length	25mm to 30mm (diff 2.0mm)



TREVALLI

BONE SCREW



Washer For Cancellous Screw

Material	Stainless Steel 316L,316LVM Titanium
Size	4.0mm, 7.0mm, 13mm
Length	Standard

DHS Lag Screw Pediatrics

Material	Stainless Steel 316L,316LVM Titanium
Size	
Length	40mm to 105mm (diff 5.0mm)



TREVALLI

INTERLOCKING NAIL SYS




K Nail for Femur

Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	07mm,08mm,09mm,10mm,11mm,12mm
Nail Length	20cm to 46cm (diff 2.0cm)




V Nail for Tibia

Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	07mm,08mm,09mm,10mm,11mm
Nail Length	20cm to 40cm (diff 2.0cm)



Rush Nail

Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	1.5mm to 6mm
Nail Length	18cm to cm (diff 2.0cm)



Supra Condylar Nail

Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	09mm,10mm,11mm,12mm
Nail Length	15cm to 30cm (diff 5.0cm)



Square Nail For Radius

Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	2mm,2.5mm,3mm,3.5mm,4mm,4.5mm,5mm
Nail Length	18cm to 30cm (diff 2.0cm)



Square Nail For Ulna

Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	2mm,2.5mm,3mm,3.5mm,4mm,4.5mm,5mm
Nail Length	18cm to 30cm (diff 2.0cm)



INTERLOCKING NAIL SYS

TREVALLi



Interlocking Tibia Nail (Solid)

Locking Bolt	4.9mm
Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	08mm,09mm,10mm,11mm
Nail Length	26cm to 40cm (diff 2.0cm)
Proximal Bolt	Locking bolt 4.9mm
Bolt Length	20mm to 60mm (diff 2.0mm)
Distal Bolt	Locking bolt 4.9mm
Bolt Length	65mm to 80mm (diff 5.0mm)

Humerus Interlocking Nail

Locking Bolt	3.1mm
Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	06mm,07mm,08mm
Nail Length	18cm to 30cm (diff 2.0cm)
Proximal Bolt	Locking bolt 3.1mm
Bolt Length	10mm to 40mm (diff 2.0mm)
Distal Bolt	Locking bolt 3.1mm
Bolt Length	10mm to 40mm (diff 2.0mm)



Humerus Interlocking Nail

Locking Bolt	4.9mm
Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	08mm,09mm,10mm,11mm,12mm
Nail Length	34.0cm to 44.0cm (diff 2.0cm)
Proximal Bolt	Locking bolt 4.9mm
Bolt Length	20mm to 60mm (diff 2.0mm)
Distal Bolt	Locking bolt 3.1mm
Bolt Length	60mm to 80mm (diff 5.0mm)

Interlocking Tibia Nail (Solid)

Locking Bolt	4.9mm
Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	08mm,09mm,10mm,11mm
Nail Length	26cm to 40cm (diff 2.0cm)
Proximal Bolt	Locking bolt 4.9mm
Bolt Length	20mm to 60mm (diff 2.0mm)
Distal Bolt	Locking bolt 4.9mm
Bolt Length	65mm to 80mm (diff 5.0mm)



INTERLOCKING NAIL SYS

TREVALLi



PFN Nail Long

Locking Bolt	3.0mm & 1.5mm
Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	9.0mm,10mm,11mm,12mm
Nail Length	25cm
Degree	130°/135°
Proximal Bolt	8.0mm
Bolt Length	40mm to 120mm (diff 5.0mm)
Proximal Bolt	6.4mm
Bolt Length	35mm to 120mm (diff 5.0mm)
Distal Bolt	Interlocking Bolt 4.9mm
Bolt Length	20mm to 60mm (diff 2.0mm)



Elastic Nail (Flexible Nail)

Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	2mm,2.5mm,3mm,3.5mm,4mm,4.5mm
Nail Length	40cm Standard



TFN Nail Short

Locking Bolt	8.0mm,6.4mm & 4.9mm
Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	9.0mm,10mm,11mm,12mm
Nail Length	18cm
Degree	130°/135°
Proximal Bolt	8.0mm
Bolt Length	40mm to 120mm (diff 5.0mm)
Proximal Bolt	6.4mm
Bolt Length	35mm to 120mm (diff 5.0mm)
Distal Bolt	Interlocking Bolt 4.9mm
Bolt Length	20mm to 60mm (diff 2.0mm)



PFN Nail Short

Locking Bolt	5.0mm,3.7mm & 3.0mm
Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	9.0mm,10mm,11mm,12mm
Nail Length	25cm
Degree	130°/135°
Proximal Bolt	8.0mm
Bolt Length	40mm to 120mm (diff 5.0mm)
Proximal Bolt	6.4mm
Bolt Length	35mm to 120mm (diff 5.0mm)
Distal Bolt	Interlocking Bolt 4.9mm
Bolt Length	20mm to 60mm (diff 2.0mm)




INTERLOCKING NAIL SYS

TREVALLi



Sirus Nail (Left & Right)

Locking Bolt	6.5mm & 4.9mm
Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	9.0mm,10mm,11mm,12mm
Nail Length	34cm to 44cm (diff 2.0cm)
Proximal Bolt	Cannulated / Solid 6.5mm
Bolt Length	30mm to 110mm (diff 5.0mm)
Proximal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm (diff 2.0mm)
Bolt Length	65mm to 80mm (diff 5.0mm)
Standard	Left & Right



Ender's Nail

Material	Stainless Steel 316L,316LVM Titanium
Nail Diameter	3mm,3.5mm,4mm,4.5mm
Nail Length	18cm to 46cm (diff 1.0cm)



TREVALLi

INTERLOCKING BOLT SYS



Inter Locking Bolt 4.9mm

Thread Diameter	4.9mm
Drill Diameter	4.0mm
Core Diameter	3.8mm
For Passing Hole	5.1mm
Material	Stainless Steel 316L,316LVM Titanium
Size	4.9mm
Length	18mm to 50mm (diff 2.0mm)
Length	65mm to 90mm (diff 5.0mm)



PFN Bolt Cannulated 6.4mm

Thread Diameter	6.4mm
Drill Diameter	5.6mm
Core Diameter	5.5mm
For Passing Hole	5.1mm
Material	Stainless Steel 316L,316LVM Titanium
Size	6.5mm
Length	18mm to 50mm (diff 2.0mm)
Length	65mm to 90mm (diff 5.0mm)



Inter Locking Bolt 3.1mm

Thread Diameter	3.1mm
Material	Stainless Steel 316L,316LVM Titanium
Size	3.1mm
Length	18cm to 50mm (diff 2.0mm)



Inter Locking Bolt 4.5mm

Thread Diameter	4.5mm
Drill Diameter	3.2mm
Material	Stainless Steel 316L,316LVM Titanium
Size	4.5mm
Length	18mm to 90mm (diff 2.0mm)



Inter Locking Bolt 3.9mm

Thread Diameter	3.9mm
Drill Diameter	3.2mm
Core Diameter	2.7mm
For Passing Hole	4.0mm
Material	Stainless Steel 316L,316LVM Titanium
Size	3.9mm
Length	18mm to 50mm (diff 2.0mm)
Length	55mm to 80mm (diff 5.0mm)



TREVALLi

INTERLOCKING BOLT SYS



PFN Bolt Cancellous 6.4mm

Thread Diameter	6.4mm
Drill Diameter	5.7mm
Material	Stainless Steel 316L,316LVM Titanium
Size	6.4mm
Length	30mm to 120mm (diff 5.0mm)



PFN Bolt Cancellous 8.0mm

Thread Diameter	8.0mm
Material	Stainless Steel 316L,316LVM Titanium
Size	8.0mm
Length	45mm to 120mm (diff 5.0mm)



PFN Bolt Cancellous 8.0mm (Hex)

Thread Diameter	8.0mm
Material	Stainless Steel 316L,316LVM Titanium
Size	8.0mm
Length	45mm to 120mm (diff 5.0mm)



TREVALLi

PROSTHESIS SYSTEM



Austin Moore Hip Prosthesis

Material	Stainless Steel 316L,316LVM Titanium
Size	37,39,41,43,45,47,49,51,53
Length	Standard



Thompson Hip Prosthesis

Material	Stainless Steel 316L,316LVM Titanium
Size	37,39,41,43,45,47,49,51,53
Length	Standard



Austin Moore Hip Prosthesis Long Stem

Material	Stainless Steel 316L,316LVM Titanium
Size	35,37,39,41,43,45,47,49,51
Length	Standard



Shoulder Prosthesis

Material	Stainless Steel 316L,316LVM Titanium
Size	35,37,39,41,43,45,47,49,51
Length	Standard




TREVALLi



Poly Axial Screw

Material	Titanium
Size	4.5mm,5.5mm,6.5mm,7.5mm
Length	20mm to 50mm (diff 5.0mm)



Mono Axial Screw

Material	Titanium
Size	4.5mm,5.5mm,6.5mm,7.5mm
Length	20mm to 50mm (diff 5.0mm)



Poly Axial Screw (D/T)

Material	Titanium
Size	4.5mm,5.5mm,6.5mm,7.5mm
Length	20mm to 50mm (diff 5.0mm)



Mono Axial Screw (D/T)

Material	Titanium
Size	4.5mm,5.5mm,6.5mm,7.5mm
Length	20mm to 50mm (diff 5.0mm)



Poly Axial Reduction Screw

Material	Titanium
Size	4.5mm,5.5mm,6.5mm,7.5mm
Length	20mm to 50mm (diff 5.0mm)



Mono Axial Reduction Screw

Material	Titanium
Size	4.5mm,5.5mm,6.5mm,7.5mm
Length	20mm to 50mm (diff 5.0mm)



TREVALLI

SPINE SYSTEM




Poly Axial Reduction Screw (D/T)

Material	Titanium
Size	4.5mm,5.5mm,6.5mm,7.5mm
Length	20mm to 50mm (diff 5.0mm)




Mono Axial Reduction Screw (D/T)

Material	Titanium
Size	4.5mm,5.5mm,6.5mm,7.5mm
Length	20mm to 50mm (diff 5.0mm)



3.5mm Dia - Polyaxial Screw

Material	Titanium
Dia	3.5mm,4.2mm
Length	12mm to 20mm (diff 2.0mm)
Length	25mm to 40mm (diff 5.0mm)



3.5mm Dia - Monoaxial Screw

Material	Titanium
Dia	3.5mm,4.2mm,5.0mm,6.0mm
Length	20mm to 50mm (diff 2.0mm)
Length	25mm to 40mm (diff 5.0mm)



Laminar Hook Narrow Blade

Material	Titanium
Dia	3.5mm,4.2mm
Length	12mm to 20mm (diff 2.0mm)
Length	25mm to 40mm (diff 5.0mm)



Laminar Hook Angle Blade

Material	Titanium
Dia	Standard
Length	Standard



TREVALLi

SPINE SYSTEM



Pedical Hook

Material	Titanium
Dia	Standard
Length	Standard



Leminar Hook Wide Blade

Material	Titanium
Dia	Standard
Length	Standard



Staple

Material	Titanium
Dia	Standard
Length	Standard



TREVALLi

SPINE SYSTEM



3.5mm Dia - Posterior Occipital Plate With Rod

Material	Titanium
Dia	3.5mm
Length	100mm, 200mm, 250mm, 300mm




Anterior Cervical Cage (Square Type)

Material	Titanium
Dia	4mm,5mm,6mm,7mm,8mm
Length	12mm x 14mm, 12mm x 16mm



Expandable Cage

Material	Titanium
Dia	10mm,12mm,14mm,16mm,18mm,20mm
Length	20mm to 75mm (diff 5.0mm)



Connecting Rod

Material	Titanium
Dia	3mm,3.5mm,5mm,5.5mm
Length	50mm to 450mm



Posterior Occipital Screw

Material	Titanium
Dia	3.5mm
Length	06mm to 20mm (Diff 2.0mm)



Eco Cage

Material	Titanium
Dia	10mm,12mm,14mm,16mm,18mm,20mm
Length	15mm to 100mm (diff 5.0mm)



TREVALLI

SPINE SYSTEM



Plif Cage

Material	Titanium
Dia	7mm,8mm,9mm,10mm,11mm
Length	20mm,25mm,30mm



Tlif Cage(Banana Type)

Material	Titanium
Dia	7mm,8mm,9mm,10mm,11mm
Length	40,80,110



Bone Screw For ACP 4mm

Material	Titanium
Dia	4mm
Length	10mm to 24mm (diff 2.0mm)



Plif Cage(Bullet Type)

Material	Titanium
Dia	8mm,10mm,12mm
Length	20mm,25mm,30mm




Occipital Loops

Material	Titanium
Dia	4mm
Length	




Anterior Cervical Plate (Z Type)

Material	Titanium
Dia	Standard
Length	25mm,27.5mm,30mm,32.5mm,35mm,37.5mm,40mm,42.5mm,45mm,47.5mm,50mm,55mm,60mm,65mm,70mm



TREVALLi

tomorros solutions . . . today!

 083 555 8643

@ trevalli@mweb.co.za

 3 Oldham Place,
Pietermaritzburg



Postnet suite #A02,
Private Bag X8,
Cascades, 3202